

# Blood Sample Procurement for Genetic Research

For submission of DNA (note address change for submission as of June 2009)

- Please collect whole blood using an EDTA or ACD Vacutainer (lavender or yellow top)
- Please include 3-5 mls of blood will provide a sufficient yield of DNA
- Samples can be shipped using STANDARD overnight (or 2-day international)
- Please ship at ambient temperature (No Ice)
- Ship samples to:

Dr. David C. Twedt  
Department of Clinical Sciences  
[College of Veterinary Medicine and Biomedical Sciences](#)  
Colorado State University  
Fort Collins, CO 80523

Also please email that a sample has been shipped: [twedt@colostate.edu](mailto:twedt@colostate.edu)

Include the following with the blood sample (Diagnostic records helpful but not required):

Reg. Name: \_\_\_\_\_ Reg. #: \_\_\_\_\_

DOB: \_\_\_\_\_ Sex: \_\_\_\_\_ Breed: \_\_\_\_\_

Sire: \_\_\_\_\_ Dam: \_\_\_\_\_

Owner's Name: \_\_\_\_\_

Address: \_\_\_\_\_

City/State/Zip: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_

**Thank you for supporting our research**

